



No barriers.
More clarity.




The ClearSight System

A simple, noninvasive approach to monitoring key hemodynamic parameters.

- Stroke Volume (SV)
- Stroke Volume Variation (SVV)
- Cardiac Output (CO)
- Systemic Vascular Resistance (SVR)
- continuous Blood Pressure (cBP)



- Finger cuff provides up to 72 hours of valuable hemodynamic insight, but continuous monitoring on one finger is limited to 8 hours. Uninterrupted monitoring over 8 hours is possible using two cuffs on two fingers.



- **Cross-section of cuff application.** To accurately mirror arterial line output, real-time finger pressure measurement is performed 1000 times per second utilizing the volume clamp method.

The ClearSight system extends continuous hemodynamic monitoring to a broader patient population, including your moderate to high-risk surgery patients. By leveraging proven ccNexfin system technology,^{1,2,5-10} the ClearSight system provides clinicians clarity without the barriers of complexity or invasiveness.

Simple for clinicians

The ClearSight system quickly connects to the patient by wrapping an inflatable cuff around the finger.¹ The simplicity of the ClearSight system gives you noninvasive access to automatic, up-to-the-minute hemodynamic information for a broader patient population.

Leverages validated technology

Used as the standard for monitoring in space for decades,³ ccNexfin noninvasive technology (volume clamp,⁴ Physiocal⁵) has been extensively validated against gold-standard monitoring technologies. Several studies have validated the ability of noninvasive ccNexfin technology to reliably track continuous blood pressure,¹ absolute cardiac output and continuous cardiac output.⁵⁻¹⁰

Advanced hemodynamic monitoring extended to moderate-risk patients

The ClearSight system enables you to expand advanced hemodynamic monitoring to your moderate to high-risk surgery patients — including elderly or obese patients^{1,11,12} — enabling you to make more informed decisions regarding volume administration.


Flow parameters provided, such as SV and SVV, have been shown to be more dynamic, and sensitive and specific for predicting fluid responsiveness than conventional, pressure-based parameters.¹³ These parameters are central to Perioperative Goal-Directed Therapy (PGDT) protocols and, when used together, are key to guiding optimal volume administration in patients at risk of developing complications.

Noninvasive simplicity. Next-generation clarity.

Edwards Lifesciences ClearSight system is comprised of the ClearSight finger cuff and EV1000 clinical platform. The system presents SV, SVV, CO, SVR and continuous BP clearly and simply on the EV1000 monitor. Easy-to-use touch screens allow selection of the parameters most meaningful for each clinical situation. This intuitive display of fundamental information minimizes unknowns and allows you to more accurately determine your patient's fluid status.

Visual clinical support

Advanced hemodynamic parameters are presented in a visually simplified manner. Color-based indicators communicate patient status at a glance, and visual clinical support screens allow for immediate recognition and increased understanding of rapidly changing clinical situations.



• Set goals and alarms

The EV1000 clinical platform allows clinicians to set specific color target ranges; and to choose parameters, alarms and targets to align with hemodynamic optimization and PGDT protocols to meet your precise patient monitoring needs.

The noninvasive ClearSight system provides valuable hemodynamic insight to an expanded patient population, for making more informed decisions about volume administration in moderate to high-risk surgery. Delivering clarity in every moment.

Historical milestones

1967

Volume clamp method to measure blood pressure invented by Jan Peñáz

1977

TNO, a high-profile Dutch research institute, starts development of noninvasive blood pressure

2005

BMEYE was founded as a spin-off from TNO

2007

Market introduction of Nexfin system with CO-Trek to non-invasively monitor blood pressure and cardiac output

2014

Market introduction of ClearSight integrated into the EV1000 clinical platform

Know More.

With the integration of the next-generation ClearSight finger cuff into the EV1000 clinical platform, a single platform can now be utilized with both noninvasive and minimally-invasive hemodynamic monitoring options, such as the FloTrac sensor. The ClearSight system also sends an analog pressure to visualize noninvasive BP on a bedside monitor. Edwards has taken this integrated system approach to protect the accuracy of the clinical information and the corresponding quality of patient care.

The ClearSight System

1. EV1000 Monitor
2. EV1000 Pump-Unit
3. Pressure Controller
4. ClearSight Finger Cuff
5. Heart Reference Sensor



Description

Description	Model No.
ClearSight Finger Cuff Small Multi Pack	CSCS
ClearSight Finger Cuff Medium Multi Pack	CSCM
ClearSight Finger Cuff Large Multi Pack	CSCL
EV1000 Clinical Platform	EV1000NI

Connectivity via IFM out through a serial connection, HL7 through an Ethernet connection or HL7 Integration Engine.

Helping to advance the care of the acutely ill for 40 years, Edwards Lifesciences seeks to provide the valuable information you need, the moment you need it. Through continuing collaboration with you, ongoing education and our never-ending quest for advancement, our goal is to delivery clarity in every moment.

Visit www.Edwards.com/ClearSight to learn more

References

1. Eeftinck Schattenkerk D, et al. Nexfin Noninvasive Continuous Blood Pressure Validated Against Riva-Rocci/Korotkoff. *American Journal of Hypertension* 2009; 22(4):378–383.
2. Martina JR, Westerhof BE, van Goudoever J, de Beauont EM, Truijen J, Kim YS, Immink RV, Jöbsis DA, Hollmann MW, Lahpor JR, de Mol BA, van Lieshout JJ. Noninvasive continuous arterial blood pressure monitoring with Nexfin®. *Anesthesiology* 2012 May;116(5):1092–103.
3. Hughson RL, Shoemaker JK, Blaber AP, Arbeille P, Greaves DK, Pereira-Junior PP, Xu D. Cardiovascular regulation during long-duration spaceflights to the International Space Station. *J. Appl. Physiol* 2012 Mar;112(5):719–27
4. Peñáz J. Photoelectric measurement of blood pressure, volume and flow in the finger. 1973; *Dresden* 1973. p.104
5. Wesseling KH, De Wit B, Van der Hoeven GMA, van Goudoever J, Settles, JJ. Physiocal, calibrating finger vascular physiology for Finapres. *Homeostasis* 1995;36:67–82
6. Kalmar AF, Vos JJ, Weening M, Mooyaart EA, Poterman M, Struys MMRF, Scheeren, TWL. Validation of Continuous Noninvasive Arterial Blood Pressure Measurements During General Anesthesia, abstract presented at the Anesthesiology Annual Meeting, 2012
7. Sterr J, Scholz S, Habicher M, Krämer M, Treskatsch S, Sander M. Comparison of the continuous noninvasive Nexfin monitoring system with conventional invasive methods to measure arterial blood pressure in high risk hip surgery, abstract presented at Euroanesthesia, 2013
8. Broch O, Renner J, Gruenewald M, Meybohm P, Schottler J, Caliebe A, Steinfath M, Malbrain M, Bein B. A comparison of the Nexfin® and transcardiopulmonary thermodilution to estimate cardiac output during coronary artery surgery. *Anaesthesia* 2012 Apr;67(4):377–83
9. Chen G, Meng L, Alexander B, Tran NP, Kain ZN, Cannesson M. Comparison of noninvasive cardiac output measurements using the Nexfin monitoring device and the esophageal Doppler. *Journal of Clinical Anesthesia* (2012) Jun 1;24(4), 275–83
10. Bubenek-Turconi SI, Craciun M, Miclea I, Perel A. Noninvasive Continuous Cardiac Output by the Nexfin Before and After Preload-Modifying Maneuvers: A Comparison with Intermittent Thermodilution Cardiac Output. *Anesth Analg*. 2013 Jun 11
11. Truijen, J et al. Noninvasive Continuous Hemodynamic Monitoring. *Journal of Clinical Monitoring and Computing* 2012;26(4):267–268
12. Maguire S, et al. Respiratory Variation in Pulse Pressure and Plethysmographic Waveforms: Intraoperative Applicability in a North American Academic Center. *Anesthesia & Analgesia* 2011;112:94–6
13. Michard F, Biais M. Rational fluid management: dissecting facts from fiction. *Br J Anaesth* 2012

For professional use. CAUTION: Federal (United States) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

Edwards Lifesciences devices placed on the European market, meet the essential requirements referred to in Article 3 of the Medical Device Directive 93/42/EEC, and bear the CE marking of conformity.

Edwards, Edwards Lifesciences, the stylized E logo, ccNexfin, ClearSight, EV1000, FloTrac and Nexfin are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

© 2014 Edwards Lifesciences Corporation. All rights reserved. AR11578

Edwards Lifesciences | edwards.com

One Edwards Way | Irvine, California 92614 USA

Switzerland | Japan | China | Brazil | Australia | India



Edwards